THUASNE®

Ortel C3

fr	Collier cervical rigide	6
en	Rigid cervical collar	7
de	Feste Halskrause	8
nl	Stijve nekbrace	10
it	Collare cervicale rigido	11
es	Collarín cervical rígido	12
pt	Colar cervical rígido	14
da	Stiv nakkekrave	15
fi	Jäykkä kaulatuki	16
sv	Rigid nackkrage	18
el	Αυχενικό κολλάρο σκληρό	19
cs	Pevný krční límec	20
pl	Kołnierz sztywny do unieruchomienia szyi	21
lv	Stingra kakla ortoze	23
lt	Standus kaklo įtvaras	24
et	Jäik emakakaela krae	25
sl	Opornica za vrat, toga	27
sk	Pevný krčný golier	28
hu	Kemény nyakmerevítő	29
bg	Твърда цервикална яка	31
ro	Guler cervical rigid	32
ru	Жесткий шейный воротник	33
hr	Tvrda ortoza za vratnu kralježnicu	35
zh	硬质颈圈	36
ar	طوق للعنق صلب	39

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegóły użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu - 使用指南 - إرشادات الاستعمال

COLLIER CERVICAL RIGIDE

Description/Destination

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Description/Destination

Anatomically-shaped chin support. ¹

Precise self-fastening height adjustment (9 cm to 13 cm). ²

Precise neck self-fastening adjustment. ³

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

Polypropylene - polyethylene - polyurethane - polyamide - polyester.

Properties/Mode of action

Stabilisation of the cervical spine.

Immobilisation of the cervical spine in severe injuries.

Indications

Cervical radiculopathy.

Cervical spondylosis.

Stable cervical fractures.

Post-operative use.

Whiplash.

Contraindications

Do not apply the product in direct contact with broken skin.

Do not use in the event of known allergy to any of the components.

Do not use in the event of traumatic injuries: severe sprain, unstable fracture.

Do not use for prolonged wear without medical supervision.

Do not use the product if the diagnosis has not been confirmed.

Precautions

Strictly comply with your healthcare professional's prescription and recommendations for use.

In the event of discomfort, significant hindrance, pain, abnormal sensations, consult a healthcare professional.

Choose the appropriate size to fit the patient, referring to the size chart.

Verify the product's integrity before every use.

It is recommended that a healthcare professional supervises the first application.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

For hygiene and performance reasons, do not re-use the product for another patient.

Do not use the device if it is damaged.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

It is recommended to adequately tighten the device to achieve support without excessive compression.

Rigid cervical collars should be removed as soon as clearance of the cervical spine is achieved, even in intubated patients or obtunded patients (with reduced level of consciousness).

Do not wear the device in a medical imaging machine.

Do not wear the device when driving a vehicle.

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

Wearing a cervical collar may cause some weakness of the neck muscles.

Wearing a cervical collar may raise intracranial pressure.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Open the self-fastener.

Place the collar around the neck. ④

The chin should rest on the chin support. ⑥

Close the self-fastener at the nape of the neck. ③

Be careful not to close the device too tightly. ⑤

Adjust the height if necessary.

Care/Maintenance

See product label. Close the self-fastening tabs before washing. Hand wash.

Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Do not tumble-dry. Do not iron. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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©Thuasne - 2047201 (2023-09)
THUASNE UK Ltd
Thuasne SAS - SIREN/RCS Nanterre 542 091 186