THUASNE®

Dorso Rigid 35

fr	Orthèse TLSO Hyperextension.....	4
en	Hyperextension TLSO braces.....	4
de	Hyperextensions-TLSO.....	5
nl	TLSO hyperextensie-orthese.....	5
it	Ortesi TLSO iperestensione.....	6
es	Órtesis TLSO hiperextensión.....	6
pt	Ortótise TLSO de hiperextensão.....	7
da	Ortose TLSO hyper-ekstension.....	8
fi	TLSO-hyperekstensiotuki.....	8
sv	Ortos för hyperextension av bröst- och ländryggen.....	9
el	Ορθωση υπερέκτασης TLSO.....	9
cs	Ortēza TLSO hyperextenze.....	10
pl	Orteza TLSO z hiper przeprostem.....	10
lv	TLSO hiperekstensijas ortoze.....	11
lt	TLSO ištempiamasis įtvaras.....	11
et	Ortoos TLSO hyper-extension.....	12
sl	Opornica TLSO za hipernategnitev.....	12
sk	Hyperextenčná TLSO ortéza.....	13
hu	TLSO ortézis hiperextenzió.....	13
bg	Лумбосакрална ортеза (TLSO) за хипер екстензия.....	14
ro	Orteză TLSO cu hiperextensie.....	15
ru	Гиперэкстензионный ортез TLSO.....	15
hr	Ortoza TLSO za hiperekstenziju.....	16
zh	TLSO超伸支架.....	16
ar	تجهيزه TLSO للمد الفائق.....	17

	cm	inch	cm	inch
S	60 - 75	24 - 29	39 - 46	15 - 18
M	75 - 90	29 - 35	42 - 49	16 - 19
L	90 - 105	35 - 41	45 - 52	18 - 20
XL	105 - 115	41 - 45	49 - 56	19 - 22

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegółowy użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu - 使用指南 - إرشادات الاستعمال

© Thuasne - 6175002 - 2025-09
© Studio Caterin

Importer in EU :
THUASNE
120, rue Marius Auphan
92300 Levallois-Perret
France
www.thuasne.com

Importer in Switzerland :
THUASNE Switzerland SA
Rue du marché 20
1204 Genève
Switzerland

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HYPEREXTENSION TLSO BRACES

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table. Available in 4 sizes.

The device is composed of:

- a frame ㉔ ensuring the brace's rigidity,
- a proximal sternal pad ㉔,
- a posterior panel ㉔ that can be positioned horizontally or vertically on the patient,
- padded lateral uprights (pads) ㉔,
- articulated and padded pelvic bar ㉔.

Composition

Rigid components: aluminum (70%) - steel (6%) - stainless steel (1%) - polyamide (8%) - polypropylene (12%) - polyoxymethylene (3%). Textile components: elastane (35%) - polyurethane (35%) - polyamide (20%) - polypropylene (10%).

Properties/Mode of action

The brace is designed to immobilise and support the thoracolumbar spine. The brace is composed of a rigid, lightweight frame comprising a proximal sternal pad, lateral uprights and a distal, articulated, padded pelvic support. It promotes spinal extension and limits forward flexion and rotation. TLSO = Thoraco-lumbo-sacral orthosis.

Indications

Postoperative immobilisation after spinal surgery. Post-traumatic - anterior body fractures (T6-L2).

Dorsolumbar / thoracolumbar osteoarthritis. Conservative treatment of stable spinal fractures (T6-L2). Osteoporotic or metastatic-related compression fractures of T6-L2. Contraindications Do not use the product if the diagnosis has not been confirmed. Do not use in the event of known allergy to any of the components. Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing. Do not use in cases of severe varicose veins. History of venous or lymphatic disorders. Preaesthesia.

Precautions

Verify the product's integrity before every use. Do not use the device if it is damaged. Choose the appropriate size to fit the patient, referring to the size chart. It is recommended that a healthcare professional supervises the first application. Strictly comply with your healthcare professional's prescription and recommendations for use. This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only. For hygiene, security and performance reasons, do not re-use the product for another patient. Wear the device above thin clothing. Ensure correct position and avoid folds when positioning the device. It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation. In the event of discomfort, significant discomfort, pain, changes in limb volume, abnormal sensations or color changes in the extremities, remove the device and consult a healthcare professional. In the event of a change in the device's performance, consult a healthcare professional. Do not wear the device in a medical imaging machine. Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Undesirable side-effect This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity. Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

HEALTHCARE PROFESSIONALS INSTRUCTION GUIDE:

Please read these instructions with the patient when fitting the brace for the first time and make sure that he/she understands how to position the brace.

Step 1: adjusting and positioning the brace

- For fitting, the patient can be supine or seated.
- Measure the patient's waist circumference at the widest point, and the height from the pelvis to just below the sternal notch (two fingers width below sternal notch) to determine the correct brace size. ①

Adjustments:

- Adjusting brace height - screw (Phillips head screwdriver).
 - Remove the screws on the lateral uprights ② and make the necessary height adjustments to adapt the brace to the patient. ①
 - Ensure that the screw on the lever latch is correctly adjusted.
 - Make sure both side uprights are adjusted to the same height.
 - Tighten all screws at the desired height.
- Adjusting brace width - screw (Phillips head screwdriver).
 - Remove the screws from the proximal sternal bar ④ and make the necessary width adjustments to adapt the brace to the patient. ①
 - Tighten all screws at the desired width.
- Posterior strap adjustment:
 - Remove the hook and loop fastener. ③
 - Adjust the length and reposition the strap. ⑦
 - Make sure the strap is securely fastened.

Fitting instructions:

Conforming the frame (especially the pelvic bar) can be done on the edge of a table or by hand. ①

- Detach frame clip ㉔ (pull the buckle string ㉔ forward to release ㉔) and open the lever latch ㉔ ㉔ ㉔.
- Position the brace on the anterior aspect of the torso. ①
- Position the proximal part of the frame two fingers width below the sternal notch.
- Wrap the attachment strap around the back, ensuring that the posterior panel is centered. ㉔
- Attach the clip to the frame. ㉔
- Once the brace is in the proper position, close the lever latch to tighten the brace. ㉔

If the brace is correctly positioned and fitted, it should remain in position and fixed despite the patient's movements.

- Remove the brace by opening the lever latch and detaching the clip from the frame.

Step 2: checking fit and alignment

- Make sure the patient can sit down and stand up without discomfort.
- If there's a gap or if the brace is too tight, adjust the brace frame or retake the pelvic circumference measurement and select the appropriate size (see size table).

PATIENT INSTRUCTIONS GUIDE:

The healthcare professional who performed the initial fitting of the brace should have made all the necessary adjustments. The healthcare professional should have also explained how to put on your brace.

- Step 1:** open the lever handle (right side) ➊ ➋ and detach the clip from the frame (left side) ➌ ➍
- Step 2:** lie down or sit down for placement.
- Step 3:** position the brace on the front of the torso. ➎
- Make sure the sternal bar is positioned approximately 2 fingers width below the sternal notch.
 - The pelvic bar is against the pelvis and positioned so that the patient can sit comfortably.
- Step 4:** wrap the strap around the back, ensuring that the posterior panel is centered. ➏
- Step 5:** attach the clip to the frame and close the lever latch to tighten the brace.
- Remove the brace by opening the lever latch and detaching the clip from the frame. ➐

Care
Product can be washed in accordance with the conditions shown on this instruction leaflet and on the label. Hand wash. Wash the rigid part with a moist cloth. Wipe the pads after each use to remove any moisture and let the pads air dry. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Dry away from any direct heat source (radiator, sun, etc.). Do not expose the product to extreme temperatures. If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Storage
Store at room temperature, preferably in the original packaging.

Disposal
Dispose of in accordance with local regulations.

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MDSS - UK RP Ltd