
THUASNE

Ortel

| | | |
|-----------|-------------------------|----|
| fr | Bandage herniaire | 4 |
| en | Hernia bandage | 5 |
| de | Bruchband | 5 |
| nl | Breukverband | 6 |
| it | Cinto erniario | 7 |
| es | Vendaje herniario | 8 |
| pt | Cinta para hérnia | 9 |
| da | Brokbind | 10 |
| fi | Tyräside | 11 |
| sv | Bräckbandage | 12 |
| el | Κηλεπίδεσμος | 13 |
| cs | Kýlní pás | 14 |
| pl | Pas przepuklinowy | 15 |
| lv | Trūces bandāža | 16 |
| lt | Išvaržos tvarstis | 17 |
| et | Songavöö | 18 |
| sl | Pas za hernijo | 19 |
| sk | Prietržový pás | 20 |
| hu | Sérvkötő | 21 |
| bg | Лента за херния | 22 |
| ro | Bandaj herniar | 23 |
| ru | Грыжевой бандаж | 24 |
| hr | Povez za herniju | 25 |
| zh | 疝帶 | 26 |
| ar | حزام الفتق | 27 |

| | | |
|--|---|--------------|
|  | 1 | 75 - 85 cm |
| | 2 | 86 - 97 cm |
| | 3 | 98 - 110 cm |
| | 4 | 111 - 125 cm |

| | |
|-----------|----------------------|
| fr | Tour de taille |
| en | Waist circumference |
| de | Tailleumfang |
| nl | Tailleomvang |
| it | Circonferenza vita |
| es | Contorno de cintura |
| pt | Perímetro de cintura |
| da | Taljemål |
| fi | Vyötärönympäry |
| sv | Midjemått |
| el | Περίμετρος μέσης |
| cs | Obvod pasu |
| pl | Obwód w talii |
| lv | Vidukla apkārtmērs |
| lt | Juostmens apimtis |
| et | Vööübermõõt |
| sl | Obseg pasu |
| sk | Obvod pása |
| hu | Derékbőség |
| bg | Обиколка на талията |
| ro | Circonferența taliei |
| ru | Обхват талии |
| hr | Opseg struka |
| zh | 腰围 |
| ar | محيط الخصر |

| | |
|-----------|----------------|
| fr | Maintien |
| en | Support |
| de | Halt |
| nl | Fixatie |
| it | Sostegno |
| es | Sujeción |
| pt | Manutenção |
| da | Støtte |
| fi | Tuki |
| sv | Stöd |
| el | Συγκράτησις |
| cs | Stabilizace |
| pl | Podtrzymywanie |
| lv | Atbalsts |
| lt | Palaikymas |
| et | Hooldus |
| sl | Oprijem |
| sk | Podpora |
| hu | Támasz |
| bg | Поддръжка |
| ro | Susținere |
| ru | Стабилизация |
| hr | Podpora |
| zh | 穩定 |
| ar | الدعم |

BANDAGE HERNIAIRE

Description/Destination

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HERNIA BANDAGE

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

Strong elastic bandage.

Removable and repositionable pads. ①

Removable and adjustable under-thigh straps. ②

Self-fastening closure system. ③

Textile components: Viscose - Cotton - Polyester - Rubber (Latex) - Elastane - Polyamid.

Non-textile components: Polyoxymethylene - Thermoplastic elastomer - Polyurethane.

Properties/Mode of action

The device ensures a support adapted to the inguinal area thanks to the removable and repositionable pads.

Indications

Bilateral ① or unilateral ② inguinal hernia.

Pre-operative wear.

Contraindications

Do not apply the product in direct contact with broken skin.

Do not use in the event of known allergy to any of the components.

Do not use on pregnant women.

Do not use in the event of non-reducible hernia.

Precautions

Verify the product's integrity before every use.

Do not use the product if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional's prescription and recommendations for use.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Do not wear the product while sleeping.

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

It is recommended to adequately tighten the device to achieve support without restricting blood circulation.

In the event of discomfort, significant hindrance, pain, abnormal sensations, remove the device and consult a healthcare professional.

For hygiene and performance reasons, do not re-use the product for another patient.

The device should not be worn long term as an alternative to surgery, as it can lead to serious complications including incarceration (trapping) of the hernia.

Contains natural latex.

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

Atrophy of the spermatic cord or fusion of the hernia sac.

Atrophy or deterioration of the fascial margins.

Strangulation of the hernia.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Preparing the device

The hernia bandage is supplied in the "bilateral hernia" position. ①

In the event of a unilateral hernia, the hernia bandage may be used for a right-side hernia or a left-side hernia.

Right-side hernia: keep the under-thigh straps and remove the left pad.

Left-side hernia: keep the under-thigh straps and remove the right pad.

Place the pad on the inside surface of the bandage against the hernia. ④

Fitting the device

Reduce the hernia before fitting the hernia bandage. To do so, it is recommended to lie on the back and massage the hernia.

Adjust the bandage to the hip circumference by positioning the fastener. ⑥

The bandage should be firmly fitted, without being too tight.

Position and adjust the length of the under-thigh straps. ⑤

Check that the pad is properly positioned against the hernia while standing.

If needed, adjust the tightening.

Care/Maintenance

Remove the pads inserts before washing. If possible use a washing net. Close the self-fastening tabs before washing. Machine washable at 40 °C (delicate programme). The device has been tested over 20 wash cycles (average: one wash per month). Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Do not tumble-dry. Do not iron. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage
Store at room temperature, preferably in the original packaging.

Disposal
Dispose of in accordance with local regulations.

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Helyezze el és állítsa be a compbántok hosszát. ©
THUASNE UK Ltd**
Thuasne SAS - SIREN/RCS Nanterre 542 091 186