
THUASNE

# LombaStart

fr	Ceinture de soutien lombaire avec compression ajustable .....	5
en	Lumbar support with adjustable compression.....	6
de	LWS-Stabilisierungsthorese mit anpassbarer Kompression.....	7
nl	Lumbale rugbrace voor de onderrug met instelbare compressie.....	8
it	Cintura di sostegno lombare a compressione regolabile .....	9
es	Faja de sujeción lumbar con compresión ajustable.....	10
pt	Cinta de suporte lumbar com compressão ajustável .....	11
da	Rygstøttebælte med justerbar kompression.....	12
fi	Lannetukivyo - säädettävä puristus .....	13
sv	Ryggördel för ryggen med justerbar kompression .....	14
el	Zώνη οσφυϊκής υποστήριξης με ρυθμιζόμενη συμπίεση .....	15
cs	Bederní podpůrný pásek s nastavitelnou kompresí .....	17
pl	Podtrzymujący pas lędźwiowy z regulowanym naciiskiem .....	18
lv	Muguras jostas atbalsta josta ar pielāgojamu kompresiju .....	19
lt	Palaikomasis reguliuojamo suspaudimo juosmens diržas .....	20
et	Nimme piirkonna reguleeritava kompressiooniga tugikorsett .....	21
sl	Oporni križni pas z nastavljivo kompresijo .....	22
sk	Podporný driekový pásek s nastaviteľnou kompresiou .....	23
hu	Ágyéki fűző állítható kompresszióval .....	24
bg	Колан за поддържане на кръста с регулируема компресия .....	25
ro	Centură de susținere lombară cu compresie reglabilă .....	26
ru	Опорный поясничный корсет с регулируемой компрессией .....	27
hr	Lumbalni potporni pojaz s podesivom kompresijom .....	29
zh	可调施压的加强型腰部稳定护带 .....	30
ar	حزام لدعم أسفل الظهر بميزة الضغط القابل للتعديل .....	31

	1	60 - 80 cm
	2	75 - 95 cm
	3	90 - 110 cm
	4	105 - 125 cm
	5	120 - 140 cm






<b>fr</b>	Tour de taille	Hauteurs : 21 cm - 26 cm
<b>en</b>	Waist circumference	Heights: 21 cm - 26 cm
<b>de</b>	Taillenumfang	Höhen: 21 cm - 26 cm
<b>nl</b>	Tailleomvang	Hoogtes: 21 cm - 26 cm
<b>it</b>	Circonferenza vita	Altezze: 21 cm - 26 cm
<b>es</b>	Contorno de cintura	Alturas: 21 cm - 26 cm
<b>pt</b>	Perímetro de cintura	Alturas: 21cm - 26cm
<b>da</b>	Taljemål	Højde: 21cm - 26cm
<b>fi</b>	Vyötärönympärys	Korkeudet: 21 cm - 26 cm
<b>sv</b>	Midjemått	Höjd: 21cm - 26cm
<b>el</b>	Περίμετρος μέσης	Υψη: 21 εκ. - 26 εκ.
<b>cs</b>	Obvod pasu	Výšky: 21 cm - 26 cm
<b>pl</b>	Obwód w talii	Wysokości: 21 cm - 26 cm
<b>lv</b>	Vidukļa apkārtmērs	Augstumi: 21 cm - 26 cm
<b>lt</b>	Juosmens apimtis	Aukštis: 21 cm - 26 cm
<b>et</b>	Võõrümbermõõt	Kõrgused: 21 cm - 26 cm
<b>sl</b>	Obseg pasu	Višina: 21 cm - 26 cm
<b>sk</b>	Obvod pása	Výšky: 21 cm - 26 cm
<b>hu</b>	Derékbőség	Magasság: 21 cm - 26 cm
<b>bg</b>	Обиколка на талията	Височини: 21 cm - 26 cm
<b>ro</b>	Circumferința taliei	Înălțimi: 21 cm - 26 cm
<b>ru</b>	Обхват талии	Высота: 21 см - 26 см
<b>hr</b>	Opseg struka	Visine: 21 cm - 26 cm
<b>zh</b>	腰围	高度 : 21 cm - 26 cm
<b>ar</b>	محيط الخصر	الارتفاعات: 21 سم - 26 سم

	Compression	Maintien	Education posturale
fr	Compression	Maintien	Education posturale
en	Compression	Support	Postural education
de	Kompression	Halt	Erlernen der Körperhaltung
nl	Compressie	Fixatie	Voorlichting over houding
it	Compressione	Sostegno	Educazione posturale
es	Compresión	Sujeción	Educación postural
pt	Compressão	Manutenção	Educação postural
da	Kompression	Støtte	Postural genopræning
fi	Kompressio	Tuki	Asentokuntoutus
sv	Kompression	Stöd	Postural träning
el	Συμπίεση	Συγκράτηση	Σωστή στάση σώματος
cs	Komprese	Stabilizace	Správné držení těla
pl	Kompresja	Podtrzymywanie	Edukacja dotycząca postawy
lv	Kompresija	Atbalsts	Stājas korekcija
lt	Kompresija	Palaikymas	Laikysenos taisymas
et	Kompressioon	Hooldus	Rühitreening
sl	Kompresija	Oprijem	Učenje telesne drže
sk	Kompresia	Podpora	Cvičenie posturálnych svalov
hu	Szorítás	Támasz	Tartásjavító oktatás
bg	Компресия	Поддръжка	Приучаване към правилна стойка
ro	Compresie	Sustinere	Educație posturală
ru	Компрессия	Стабилизация	Исправление осанки
hr	Kompresija	Potpore	Učenje o držanju
zh	压紧	稳定	姿势联系
ar	الضغط	الدعم	تدريب وضعی

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en

## **LUMBAR SUPPORT WITH ADJUSTABLE COMPRESSION**

### **Description/Destination**

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Elastic fabric. ①

Abdominal panels with self-fastening surface, finger-loop and thumb-loop. ②

Doubled additional strap with finger-loop and thumb-loop. ③

Conformable dorsal stays. ④

### **Composition**

Textile components: polyamide - polyester - elastane - polyurethane.

Non-textile components: brushed tempered steel - polyoxymethylene - polypropylene - polyethylene.

### **Properties/Mode of action**

The device relieves low back pain by supporting the spine through an increased intra-abdominal pressure which eases the load on the intervertebral discs.

### **Indications**

Non-specific low back pain (acute, subacute and chronic phases). Prevention of recurrent low back pain in combination with return to activity and/or therapeutic information.

Spinal stenosis.

Spine stability.

Resumption of activity.

### **Contraindications**

Do not use in the event of known allergy to any of the components.

\*Do not use on pregnant women.

Do not use in the event of a hiatus hernia.

Do not use in the event of bone cancer with spinal metastases.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom an increase in arterial blood pressure is not recommended.

### **Precautions**

Verify the product's integrity before every use.

Do not use the product if it is damaged.

Wear the device above thin clothing.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

Do not wear the product while sleeping.

In the event of discomfort, significant hindrance, pain, abnormal sensations, remove the device and consult a healthcare professional.

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

Measure the waist circumference of the patient at the mid-distance between the iliac crests and floating ribs.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a healthcare professional supervises the first application.

For hygiene, security and performance reasons, do not re-use the product for another patient.

Do not wear the product in a medical imaging machine.

### **Undesirable side-effect**

This device can cause skin reactions (redness, itching, irritations...).

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

### **Instructions for use/Application**

It is recommended that a healthcare professional adjust the stays to the shape of the patient's back before the belt is worn for the first time. ⑤

Position the belt around the waist, using abdominal panel and the finger loop located on the opposite abdominal panel. ⑥

The Thuisne logo indicates which way up the belt should be put on.

Pull the finger loop and the opposite abdominal panel simultaneously and evenly.

Fasten the belt, right side over left side, adjusting to the desired compression. ⑦

The belt must be centred on the spine: the dorsal stays should be positioned on either side of the spine. ⑧

The bottom of the belt must stop at the top of the crease between the buttocks. ⑨

For a reinforced compression, use the additional strap. Pull on the additional strap simultaneously and evenly. ⑩

Attach the self-fastening parts of the additional strap onto the belt's abdominal panels or place them on top of each other. ⑪

The level of compression can be adjusted as required throughout the day.

## Care/Maintenance

Product can be washed in accordance with the instructions shown on this leaflet and on the label. Close the self-fastening tabs before washing. Machine washable at 30°C (delicate programme). It is not recommended to machine wash the device more than 15 times (average: one wash per month). Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Do not tumble-dry. Do not iron. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

## Storage

Store at room temperature, preferably in the original packaging.

## Disposal

Dispose of in accordance with local regulations.

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La cintura deve essere centrata sulla colonna vertebrale: le stecche dorsali devono trovarsi su entrambi i lati della colonna vertebrale. ©
**THUASNE UK Ltd**
Thuasne SAS – SIREN/RCS Nanterre 542 091 186