
THUASNE®

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

fr	Genouillère pour Arthrose	4
en	Knee Osteoarthritis brace	4
de	Orthese für Kniearthrose	5
nl	Kniebrace voor artrose in de knie	6
it	Ginocchiera per gonartrosi	7
es	Rodillera para Artrosis de rodilla	8
pt	Joelheira para Artrose do joelho	8
da	Knæskinne til artrose i knæet	9
fi	Polviartroosin polvituki	10
sv	Knästöd för knärartros	11
el	Επιγονατίδα για οστεοαρθρίτιδα του γόνατος	11
cs	Kolenní ortéza při artróze kolene	12
pl	Orteza stawu kolanowego przy chorobie zwyrodnieniowej stawu kolanowego	13
lv	Ceļa ortoze ceļgala osteoartrītam	14
lt	Kelio įtvartas kelio osteoartritui	15
et	Põlvetugi põlveliigese artroosi korral kasutamiseks	15
sl	Kolenska opornica za artrozo kolena	16
sk	Kolenná ortéza na osteoartrózu kolena	17
hu	Térdrögzítő térdízületi artrózis esetén	18
bg	наколенка при остеоартрит на коляното	18
ro	Genunchieră pentru artroza genunchiului	19
ru	Наколенник для лечения артроза колена	20
hr	Ortoza za koljeno kod artroze koljena	21
zh	膝骨关节炎用护膝	22
ar	دعامة ركبة خاصة بفصال الركبة	23


A


OU - OR - ODER - OF - O - O - OU - ELLER - TAI - ELLER - 'H - NEBO - LUB - VAI - ARBA - VÖI - ALI - ALEBO -
VAGY - ИЛИ - SAU - ИЛИ - ILI - 或 - 9

B

	1 A	41 - 45 cm
	2 A	45 - 50 cm
	3 A	50 - 55 cm
	4 A	55 - 61 cm

	1 B	41 - 45 cm
	2 B	45 - 50 cm
	3 B	50 - 55 cm
	4 B	55 - 61 cm

fr	Circonférence de la cuisse (15 cm au-dessus du genou)
en	Thigh circumference (15 cm above the knee)
de	Oberschenkelumfang (15 cm über dem Knie)
nl	Omvang van het dijbeen (15 cm boven de knie)
it	Circonferenza della coscia (15 cm al di sopra del ginocchio)
es	Circonferecia del muslo (15 cm por encima de la rodilla)
pt	Circunferência da coxa (15 cm acima do joelho)
da	Lårets omkreds (15 cm over knæet)
fi	Reiden ympärys (15 cm polven yläpuolelta)
sv	Omkrets runt lår (15 cm ovanför knäet)
el	Περιμετρος του μηρού (15 εκ. πάνω από το γόνατο)
cs	Obvod stehna (15 cm nad kolenem)
pl	Obwód uda (15 cm powyżej kolana)
lv	Augšstilba apkārtmērs (15 cm virs ceļgala)
lt	Šlaunies apimtis (15 cm virš kelio)
et	Reie ümbermõõt (15 cm põlvest kõrgemal)
sl	Obseg stegna (15 cm nad kolenom)
sk	Obvod stehna (15 cm nad kolenom)
hu	Combókőrfogat (15 cm a térd felett)
bg	Обиколка на бедрото (15 cm над коляното)
ro	Circumferința coapsei (la 15 cm deasupra genunchiului)
ru	Диаметр окружности бедра (на 15 см выше колена)
hr	Opseg bedra (15 cm iznad koljena)
zh	大腿周长 (膝盖上方15厘米处)
ar	محيط الفخذ (15 سم فوق الركبة)

*TensionStrap System

fr	Compression
en	Compression
de	Kompression
nl	Compressie
it	Compressione
es	Compresión
pt	Compressão
da	Kompression
fi	Kompressio
sv	Kompression
el	Συμπίεση
cs	Komprese
pl	Kompresja
lv	Kompresija
lt	Kompresija
et	Kompressioon
sl	Kompresija
sk	Kompresia
hu	Szorítás
bg	Компресия
ro	Compresie
ru	Компрессия
hr	Kompresija
zh	压紧
ar	الضغط

fr	Décharge - *système de force de rappel élastique
en	Off-loading - *elastic resistance force system
de	Entlastung - *elastisches Rückstellkraftsystem
nl	Losmaken - *elastisch retourkrachtsysteem
it	Scarico - *sistema di forza di ritorno elastica
es	Descarga - *sistema de fuerza de retroceso elástico
pt	Alívio do peso - *sistema de força elástica de retorno
da	Aflastning - *elastisk system til styrkegenoprettelse
fi	Vapautus - *elastinen palautusvoimajärjestelmä
sv	Avlastning - *elastiskt korrigeringsssystem
el	Αποφόρτιση - *σύστημα ελαστικής δύναμης επαναφοράς
cs	Snížení zátěže - *systém elastické zpětné síly
pl	Odciążenie - *system gumek zapewniających siłę powrotną
lv	Atslogošana - *elastīgā atgriezes spēka sistēma
lt	Apkrovos sumažinimas - *tampri tamprumo jėgos sistema
et	Koormuse vähendamise - * elastne tagasitõmbesüsteem
sl	Razbremenitev - *sistem elastične povratne sile
sk	Odlahčenie - *elastický upínací systém
hu	Telhermentesítés - *rugalmas visszatérítő erő rendszer
bg	Отбременяване - *система на базата на еластична връщаща сила
ro	Descărcare - *sistem elastic bazat pe forța de revenire
ru	Разгрузка - *система эластичной силы сопротивления
hr	Rasterećenje - *elastični sustav za učvršćivanje
zh	减压 - *弹性回归力系统
ar	تخفيف ثقل - *نظام قوة الارتداد المطاطية

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Description/Destination

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KNEE OSTEOARTHRITIS BRACE

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Available in 4 sizes and 2 versions:

- Medial arthrosis on the left / lateral arthrosis on the right **A**.
- Medial arthrosis on the right / lateral arthrosis on the left **B**.

The device is composed of:

- an anatomically-shaped compressive knit.
- 3-point TensionStrap System (elastic resistance force system).
- a flexible stay.

Composition

Textile components: polyamide - polyester - elastane - polyurethane.

Zinc and steel stay.

Properties/Mode of action

X-shaped 3-point off-loading system, with a TensionStrap System (Elastic resistance force system) to ensure alignment and off-loading of the affected compartment.

An anatomically-shaped knit for even compression.

Indications

Symptomatic (mild) unicompartmental femorotibial osteoarthritis (arthrosis).

Contraindications

Do not use in the event of known allergy to any of the components.

Do not apply the product in direct contact with broken skin.

Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Precautions

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size, referring to the size chart.

For hygiene and performance reasons, do not re-use the product for another patient.

In case of history of venous or lymphatic disorders, consult a healthcare professional.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

Before any sports activity, check the compatibility of the use of this medical device with your healthcare professional.

Consult a healthcare professional in case of uncertain diagnosis.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended to wear the device directly on the skin, unless contraindicated.

- Slip the knee support on the leg with your foot extended ①.

- Make sure the kneecap is positioned in the center of the ribbed area provided ②, and that the elastic straps cross over on the side opposite the painful area.

To remove the brace: grasp the bottom of the brace to remove it.

Care

Machine washable at 30°C (delicate programme). If possible use a washing net. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Do not tumble-dry. Do not iron. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). It is not recommended to machine wash the device more than 24 times (average: two washes per month). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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©Thuasne - 2055201 (2024-08)
THUASNE UK Ltd
Thuasne SAS – SIREN/RCS Nanterre 542 091 186